

LS SCIENTIFIC®

Laboratory solutions experts



Solutions Provider of
Laboratory Products
and Services



ABOUT US

At LS Scientific, we believe in delivering more than just products or services — we deliver trust, quality, and long-term value. Founded on a commitment to excellence and customer satisfaction, we have grown into a trusted partner for clients globally.

Our team combines industry expertise with a passion for innovation, ensuring that every solution we offer is tailored to meet the unique needs of our customers. From our attentive service to our high standards, everything we do is focused on helping you succeed.

Whether you're looking for reliable products, expert advice, or a partner who truly understands your goals, LS Scientific is here for you — today and into the future.



BUSINESS GOALS & OBJECTIVES

To ensure superior standards and enhance the quality of finished products, using innovative technology to satisfy global scientific requirements.



OUR VALUES

Our customers are kings, and we are passionate about maintaining the high standard of quality we preach. As such, we infuse honesty, integrity, and trust into all aspects of our business.



OUR MISSION

To rid the world of substandard products and improve Quality of Life.



OUR PRODUCTS

Laboratory Chemicals

Our extensive chemical inventory caters to laboratories of all sizes and research needs, offering a diverse selection from analytical-grade to high-purity HPLC-grade chemicals, as well as essential reference standards.



Laboratory Equipment & Instrument

Our equipment portfolio offers a diverse range of instruments, catering to both novice and experienced users and providing solutions for a variety of analytical needs.



Glassware & Consumables

Our glassware offering features a wide range of high-quality products at competitive prices, ensuring that laboratories and research facilities have access to reliable and affordable equipment.



Safety & PPE

Our comprehensive safety equipment ranges from personal protective gear like gloves and bouffant caps to specialised containment solutions such as fume cabinets and laminar flow cabinets.



**NEW
PRODUCT!**

LS SCIENTIFIC®

BottleTop Dispenser

Precision Made Simple

Easy Re-Calibration: The calibration tool is integrated and can be easily recalibrated using the position indicator

Accurate Volume Control: Set volumes with a clear, large-scale display and a secure slide-lock mechanism

Splash Protection: No splashes during high-pressure dispensing

Multi-Bottle Fit: Adapters for varied neck sizes for different bottles

Effortless Dispensing Mechanism: Low-resistance valve ensures smooth piston movement

Ordering Information

Product code	Description	Volume range	Increment (ml)	Accuracy ± %	Accuracy ± ml	CV ± %	CV ± ml
LSBTDP0025	LS-BTD-Precision 2.5 - Bottle Top Dispenser 0.25ml - 2.5ml, with Recirculation Valve	0.25 - 2.5	0.05	0.6	0.012	0.2	0.004
LSBTDP0050	LS-BTD-Precision 5.0 - Bottle Top Dispenser 0.50 - 5.0ml, with Recirculation Valve	0.5 - 5.0	0.1	0.6	0.03	0.2	0.01
LSBTDP0100	LS-BTD-Precision 10.0 - Bottle Top Dispenser 1.0-10ml, with Recirculation Valve	1.0 - 10.0	0.2	0.6	0.06	0.2	0.02
LSBTDP0250	LS-BTD-Precision 25 - Bottle Top Dispenser, 2.5 - 25ml with Recirculation Valve	2.5 - 25.0	0.5	0.6	0.18	0.2	0.06
LSBTDP0500	LS-BTD-Precision 50 - Bottle Top Dispenser 5.0 - 50ml, with Recirculation Valve	5.0 - 50.0	1	0.6	0.36	0.2	0.12
LSBTDP1000	LS-BTD-Precision 100 - Bottle Top Dispenser 10.0 - 100ml, with Recirculation Valve	10.0 - 100.0	2	0.6	0.6	0.2	0.2



OUR >>>>>
STRATEGIC
PARTNERS



Chemicals & Solvents

Based in Belgium, Chem-Lab is a trusted supplier of high-purity chemicals and certified reference materials. As part of AnalytiChem, Chem-Lab specialises in the production of tailor-made reagents, enabling us to offer customised chemical solutions to meet your specific laboratory or industrial requirements. All products are manufactured in ISO-certified facilities and are supplied with comprehensive documentation, including certificates of analysis, ensuring full traceability and compliance.



Analytical Equipment & Instruments

Velp Scientifica is a recognised leader in the development of laboratory instruments designed to enhance productivity and precision in scientific and industrial settings. Their extensive portfolio includes stirring and mixing equipment, digestion and distillation systems, and a range of analytical solutions. Known for their innovation and reliability, Velp's technologies—such as those in the UDK Series—are built to meet demanding international standards, offering flexible automation levels to support a variety of laboratory applications and throughput needs.



Glassware & Plasticware

DWK Life Sciences is a leading manufacturer of high-quality laboratory consumables and packaging solutions, supporting scientific and industrial workflows across the globe. With trusted brands such as WHEATON®, AZLON®, PYREX®, QUICKFIT®, and MBL®, DWK offers a comprehensive range of glass and plastic products known for precision, reliability, and performance. Whether in research, quality control, or production environments, DWK provides the tools professionals depend on to drive results with confidence.



designed for scientists

Liquid Handling, Separation & Mixing



IKA is a global leader in the design and manufacture of innovative laboratory and process equipment, trusted by scientists and engineers worldwide. Their advanced portfolio includes centrifuges, viscometers, mills, pipettes and rotary evaporators —engineered to deliver precision, reliability, and ease of use. Renowned for their German craftsmanship and commitment to innovation, IKA's solutions are designed to meet the highest standards, supporting research, development, and production processes across diverse scientific and industrial applications.

Technica



Carbonation

Technica is a specialist in gas dosing technology designed for the food and beverage industry. Their innovative systems support beverage formulation and development, enabling precise control over carbonation and dosing processes. Technica's solutions are widely used to help optimise product quality and ensure consistency in sensory characteristics across a range of beverage applications.



OUR ▶▶▶▶▶ PRODUCT CATEGORIES

Autoclaves



Balance & Weighing



Cabinets



Calibrated Test Weights



Carbonation Systems



Chemicals



Colony Counters



Digestion



Dispensing



Distillation



Electrodes & Probes



Glassware



Hot Plates



Incubators & Ovens



Lab Tools



Microbiology



Microscopes



Mills



Mixers



Plasticware



Reagents



Refrigeration &
Freezers



Safety & PPE



Sample & Liquid
Handling



Separation



Shakers



Stirrers



Test Kits & Strips



Testing Equipment



Thermometers



Titration



Water Bath



Water Purification
Systems





Unit 112 Joseph Wilson Industrial Estate, Millstrood Road,
Whitstable, Kent, CT5 3SN
United Kingdom

🌐 www.lsscscientific.com

📘🌐 LS Scientific 📺 @lsscscientific1510

☎ +44 (0) 1227 200 161 ✉ uksales@lsscscientific.com